<p>REF CVRT2500 (25 tests)</p>	<h1>nCOVID-19 IgG&amp;IgM POCT</h1>		
<p><b>INSTRUCTIONS FOR USE</b> EN English</p>		  <p>25</p>	

**[Product Name]**

nCOVID-19 IgG&amp;IgM POCT

**[Package Size]**

25 tests/kit

**[Intended Use]**

The nCOVID-19 IgG&IgM POCT (point-of-care testing) kit is intended to be used as a diagnostic aid when assessing the immunological status of patients with suspected coronavirus disease (COVID-19).

The kit is used for the qualitative detection of the presence of anti-SARS-COV-2 immunoglobulin G (IgG)/ immunoglobulin M (IgM) in human serum, plasma, and whole blood, using solid phase colloidal gold immunochromatographic technology.

NOTE: This product must be used by medical professionals in strict compliance with the instructions in this document.

CAUTION: Medical decisions should not be based exclusively on the result of this test; instead, all available clinical and laboratory data must be taken into account as a whole.

**[Clinical significance]**

The 2019 novel coronavirus, recently named SARS-CoV-2 (originally called 2019-nCoV) was identified as the causative agent of an outbreak of viral pneumonia, COVID-19.

SARS-CoV-2 belongs to the  $\beta$ -coronavirus genus and is similar to the SARS in 2003 and MERS in 2012. The genome of coronavirus encodes four structural proteins including Spike (S) protein, Envelope (E) protein, Membrane (M) protein, and Nucleocapsid (N) protein. The body's immune system generally produces specific antibodies 1 to 2 weeks after the initial SARS-CoV-2 infection.

**[Principles of Procedure]**

The nCOVID-19 IgG&IgM POCT kit adopts the solid phase colloidal gold immunochromatographic technology for the qualitative determination of the presence of IgM/IgG antibodies against SARS-CoV-2 in human serum, plasma, and whole blood. The gold SARS-CoV-2 antigen conjugate and the gold chicken IgY conjugate are coated to the conjugate pad in advance. The test line T1 (antibodies against human IgM), the test line T2 (SPA) and the control line (antibodies against chicken IgY) are pre-coated on the surface of Nitrocellulose (NC) membrane.

When the specimen is added to the sample pad, it migrates through the conjugate pad. In the presence of sufficient IgM antibodies against SARS-CoV-2, a gold SARS-CoV-2 antigen conjugate - IgM antibodies against SARS-CoV-2- antibody against human IgM complex is formed, resulting in the appearance of test line T1 on the strip, i.e. IgM Positive. Further, in the presence of sufficient IgG antibodies against SARS-CoV-2, a gold SARS- CoV-2 antigen conjugate - IgG antibodies against SARS-CoV-2 – SPA complex is formed, resulting in and the appearance of test line T2 on the strip, producing an IgG positive result. If the specific IgM/IgG antibodies are absent, or present at a very low level, no test line appears, showing a Negative result.

**[Materials Provided]**

1. Test cassette: The gold SARS-CoV-2 antigen conjugate and the gold chicken IgY conjugate are coated to the conjugate pad; the antibodies against human IgM, the SPA and the antibodies against chicken IgY are pre-coated on the surface of Nitrocellulose (NC) membrane. 1×25 cassettes

2. Sample diluent	4ml x1 bottle
3. Safety lancet	1×25 pieces
4. Alcohol pad	1×25 pieces
5. Disposable transfer pipette	1×25 pieces
6. Instructions for use	1 piece

### [Storage Requirements and Shelf-life]

The test cassettes and the sample diluent must be stored between 4-30°C for up to 24 months. Do not use frozen devices. Do not use after expiry date. Refer to the packaging label for production date and expiry date.

### [Specimen types]

- (1) Human serum, plasma, and whole blood (venous and capillary) specimens are validated to be used with this assay.
- (2) Serum and plasma specimens may be stored up to 8 hours at 18-26°C or 2-8°C for up to 7 days from time of blood drawing; or be stored at -20°C or below for long time storage. Several freeze-thaw cycles should be avoided.  
  
Venous whole blood specimens containing anticoagulants can be stored up to 8 hours at 18-26°C or up to 7 days at 2-8°C and should not be frozen.  
  
Whole blood specimens collected by fingerstick should be used freshly.
- (3) Plasma specimens can be prepared with EDTA, heparin or sodium citrate as anticoagulant.

### [Specimen Collection and Storage Requirements]

#### [Whole Blood]

- Collection by venipuncture (Whole Blood with anticoagulant)

Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA or sodium citrate) by venipuncture.

Venous whole blood specimens containing anticoagulant can be stored up to 8 hours at 18-26°C or up to 7 days at 2-8°C if not tested at the same day of specimen collection and should not be frozen. The specimens should be brought to room temperature (18- 26°C) prior to use.

- Collection by fingerstick (Whole Blood without anticoagulant)


Figure a

Figure b

Figure c

Figure d

1. Choose the ring finger of the patient. Massage and/or shake to stimulate blood flow towards the collection area.
2. Open the alcohol pad and clean the fingertip (Figure a). Allow the fingertip dry.
3. Remove the cap of the safety lancet. (Figure b)
4. Place the lancet on the fingertip, press the lancet against the finger and puncture the skin. (Figure c)
5. Discard the lancet into a bio-hazard waste container immediately.
6. Squeeze the fingertip to produce a bulb on the skin's surface, gently place the disposable transfer pipette on the fingertip and collect the drop of the blood. (Figure d)
7. Allow the blood to reach the first graduation and the specimen volume will be 20µL. (Figure e)
8. Avoid formation of bubbles during the collection of specimens. Disposable transfer pipette is used for fingerstick only.

Figure e

9. Transfer the specimen to the sample area immediately after collection. Avoid formation of bubbles.
10. Discard the disposable transfer pipette into a bio-hazard waste container immediately.

### [Plasma or Serum]

- Plasma

Draw whole blood into the collection tube (containing anticoagulants such as heparin, EDTA, or sodium citrate) by the venipuncture technique and lightly agitate by hand, then centrifuge the whole blood sample to produce a plasma supernatant.

- Serum

Draw the whole blood into the collection tube by venipuncture, leave to settle for 2 hours for blood coagulation and then centrifuge blood to get serum supernatant.

Serum and plasma specimens may be stored up to 8 hours at 18-26°C or 2-8°C for up to 7 days from time of draw, or be stored at -20°C or below for long time storage. Several freeze-thaw cycles should be avoided.

Plasma or serum specimens should be adjusted to room temperature prior to use.

Plasma or serum specimens with visible bio-contamination should not be used.


Haemolytic or turbid specimens should not be used.

### [Test Procedure]

PLEASE USE IN STRICT ACCORDANCE WITH THE INSTRUCTIONS FOR USE.

Equilibrate all specimens and the devices to room temperature before testing. Sample diluent from different batches cannot be mixed.

1. Take out a test cassette from the foil pouch, and place it on a horizontal surface. Test cassette should be used immediately upon opening the pouch.
2. Add 10 µL of plasma or serum specimen or 20 µL whole blood specimen to the sample area firstly, then slowly add 3 drops of sample diluent (120 µL) to the same area. In case of fingerstick specimen, sample diluent should be added immediately after the sample is added into sample area.
3. Incubate the cassette at room temperature (18-26°C) and read the result after 15 minutes but not more than 20 minutes. Reading the test before 15 minutes or after 20 minutes may cause false result.



Add 10µL of Plasma or serum specimen

Add 20µL whole blood specimen

**[Interpretation of Results]**

T1: IgM test line; T2: IgG test line; C: control line

Positive	1 IgM Positive/IgG Negative: Visible reddish-purple bands appear both at the control line (C- line) and the test line 1 (T1) of the cassette.	C T2 T1 
	2 IgG Positive/IgM Negative: Visible reddish-purple bands appear both at the control line (C- line) and the test line 2 (T2) of the cassette	C T2 T1 
	3. IgM Positive/IgG Positive: Visible reddish-purple bands appear at the control line (C-line), both the test line 1 (T1) and the test line 2 (T2) of the cassette.	C T2 T1 
Negative	A visible reddish-purple band appears only at the control line (C-line) of the cassette	C T2 T1 
Invalid Result	No visible reddish-purple band appears at the control line or the test line of the cassette, or a reddish-purple band appears only at the test line (T-line) of the cassette	C T2 T1    C T2 T1    C T2 T1    C T2 T1 

Interpretation of the results :

- 1)When the test results show “negative”, but the patient has associated symptoms, additional testing should be carried out after a few days to ensure that no antibodies have formed.
- 2) When the test results show “positive”, further testing should be carried out to verify the result.

**[Limitations of the Procedure]**

1. The kit is designed to detect IgM/IgG antibodies against SARS-CoV-2 in human serum, plasma, and whole blood. Specimens other than the specified ones may not produce accurate results.
2. The nCOVID-19 IgG&IgM POCT kit uses solid phase colloidal gold immunochromatographic technology for the qualitative determination only
3. This product is intended for detection of antibodies against SARS-CoV-2 in individuals to aid in the diagnosis of COVID-19 infection. The final results should be evaluated by medical doctors and can be used for clinical reference but should not be the only basis for the diagnosis and treatment. The clinical management of patients should be considered in combination with patients’ symptoms and medical history, other laboratory tests, treatment response, epidemiology and other information. Negative results with the antibody tests do not preclude SARS-CoV-2 infection.
4. In case of doubtful results due to manual operations and sample collection, repeated testing should be done to ensure consistent results.
5. Cross-reactivity was found in IgM+/IgG- clone specimen, which resulted in the false positive result of IgG test line, but cross-reactivity has not been verified in genuine specimens yet. Although it suggests the potential interference of IgM to the IgG test line, it will not affect the interpretation of SARS-CoV-2 total antibodies.

**[Warnings]**

1. This product is only used for in-vitro diagnostic purposes; the manual should be read in detail before use.
2. Incorrect reporting and false results may occur if not used in strict accordance with the instructions for use.
3. Do not use the test cassette if the pouch is damaged or the seal is broken.
4. Do not use the accessories in the kit if the accessories are damaged or broken.
5. The kit assay should be considered as infectious. The infectious disease laboratory inspection procedure should be

followed.

6. Ensure accurate quantities of samples for testing; too little/too much sample size may lead to incorrect results.
7. This product requires visual reading. In order to ensure correct interpretation of the results, do not read in a dark environment.
8. Sample diluent should be added immediately after the sample is added into sample area. Delay of adding sample diluent may lead to invalid test results.

### [Performance characteristics]

#### 1. Diagnostic specificity and sensitivity:

The following diagnostic specificity and sensitivity data was obtained from a clinical study conducted on the Italian population:

- Diagnostic specificity

The specificity study was worked out on 357 samples collected from 51 patients with negative PCR (of which 3 positive for FluA and 3 for FluB), 30 asymptomatic patients of COVID-19 period, 176 blood donors before 2019 and 100 patients of TORCH laboratory routine collected in 2017 and 2018, who underwent testing via nCOVID-19 IgG&IgM POCT. The results were anti-SARS-CoV-2 IgG/IgM negative in 355/357 samples, yielding a specificity of 99.4% (C.I. at 95%: 97.8% – 99.9%).

- Diagnostic sensitivity

148 specimens of PCR positive patients were tested. Depending on the time of onset of the symptoms, the following values of antibody positivity were obtained:

From (days)	To (days)	Results		
		POS. IgG/IgM POCT (N)	POS. PCR (N)	Sensitivity (%)
0	5	4	22	18.2
6	10	22	40	55.0
11	15	25	28	89.3
16	40	58	58	100.0

Temporal changes in antibody positivity are reflected in the table.

Considering a seroconversion period of  $\geq 16$  days, 58 of 58 PCR positive results were anti-SARS-CoV-2 IgG/IgM positive via nCOVID-19 IgG&IgM POCT, yielding a sensitivity of 100.0% (C.I. at 95%: 92.3% -99.8%).

#### 2. Analytical specificity

- Endogenous interference:

8 endogenous interfering substances were tested. For each interfering substance, no obvious influence was found in plasma or in whole blood, caused by these 8 interfering substances, when they are at the concentration as below.

#### Concentrations of endogenous interfering substances

Endogenous interfering substances	Concentrations planned to be verified
Bilirubin Conjugated	5 mg/dL
Bilirubin Unconjugated	15 mg/dL
Hemoglobin	200 g/L
Human Serum Albumin	50 g/L
Intralipid	5.6 mmol/L
Cholesterol	6.5 mmol/L
Rheumatoid Factor (RF)	2000 IU/ml
Anti-nuclear Antibodies (ANA)	Positive at 1:1280 dilution

The study has shown that the performances of the kit are not influenced by the presence in the specimen of the potentially

endogenous interfering substances listed in the table above, up to the tested concentration

• Cross-reactivity evaluation:

148 potential interfering specimens including 17 types were tested to evaluate the cross reactivity. The clinical cross interference specimens include: human coronavirus 229E (5), human coronavirus OC43 (2), human coronavirus NL63 (3), MERS coronavirus (3), influenza A (FluA) (15) and B (FluB) (14), respiratory syncytial virus A (RSVA) (5) and B (RSVB) (3), Rhinovirus (RV) (6), Mycoplasma pneumoniae (MP) (10), Enterovirus (EV) (3), human metapneumovirus (HMPV) (5), Varicella zoster virus (VZV) (10), Chlamydia pneumoniae (CPN) (34), Cytomegalovirus (CMV) (10), Human immunodeficiency virus (HIV) (10), Anti-nuclear Antibodies (ANA) (10) . There are no cross reactions found in 147 specimens. Among the 34 Chlamydia pneumoniae specimens, 33 specimens were nonreactive and 1 specimen was reactive.

3. Repeatability and reproducibility

Two different operators tested 4 specimens (including IgM low reactive, IgG low reactive, IgG medium reactive and negative specimen) independently with a same batch of products in two different sites, continuously for 5 days. Each specimen was repeated 10 tests in a run.

The expected negative specimens showed negative results and the positive specimens showed positive results in all tests. The intensity of T1 and T2 lines of a same specimen was consistent and uniform and did not show any obvious difference among all measures.

• Repeatability:

The study demonstrated the repeatability of the method.



• Reproducibility:

The study demonstrated the reproducibility of the method.

**[References]**

1. National health committee “new coronavirus infection pneumonia diagnosis and treatment plan” (fifth edition)
2. World Health Organization: Clinical management of severe acute respiratory infection when Novel coronavirus (nCoV) infection is suspected: Interim Guidance

**[Key to symbols used]**

	<b>Caution</b>		<b>Temperature limitation (4-30 °C)</b>
	<b>Batch code</b>		<b>Date of manufacture</b>
	<b>In vitro Diagnostic use Medical device</b>		<b>Consult instructions for use</b>
	<b>Catalogue number or order number</b>		<b>Do not use if package is damaged</b>
	<b>Use by</b>		<b>Sufficient for “n” use</b>
	<b>Do not reuse</b>		<b>Biological risks</b>
	<b>Sterilized using irradiation</b>		<b>Sterilized using Ethylene Oxide</b>
	<b>Manufacturer</b>		<b>CE marking</b>

TECHNOGENETICS S.r.l.  
Corso Vittorio Emanuele II, 15  
20122 Milano – Italia  
[www.technogenetics.it](http://www.technogenetics.it)

Updates are shown on a gray background

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